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Dentaurum Implants GmbH I Turnstr. 31 I 75228 Ispringen I Germany

2. Brief description

One CAD/CAM Titanium base, CEREC® is supplied with one abutment screw, both in non-sterile condition. All components are intended for single use only. The scan body has to be ordered separately at Sirona. Both individually manufactured mesostructures and temporal restorations can be glued on the CAD/CAM Titanium base, CEREC®. The glued components are screwed to the implant with the abutment screw in the patient's mouth. The scan body's sole purpose is to map the implant position for creating the design in the inLab 3D software.

The CAD/CAM Titanium base, CEREC® is available in different versions, each being compatible with a specific diameter of the tiologic® implant system.

Type of implant	Implant diameter	CAD/CAM Titanium base by Dentaurum Implants	TiBase scan body by Sirona
tiologic®	3.3 mm	386-195-00	62 82 615
	3.7 mm	386-295-00	
	4.2 mm	386-395-00	
	4.8 mm		
	5.5 mm		
tioLogic® ST	3.3 mm	386-195-00	62 82 615
	3.7 mm	386-295-00	
	4.2 mm	386-395-00	
	4.8 mm		
	5.5 mm		

3. Materials

CAD/CAM Titanium base	3.7165
AnoTite screw	3.7165

4. Intended use

The Sirona Dental CAD/CAM system is indicated for optical impressions to record the topographical properties of teeth, dental impressions or plaster models by computer-aided design and manufacturing in patients that need dental restorations or prosthetics. The system can also be used to process mesostructures (restorative dental prosthetics used in combination with endosseous implant abutments). The system that can be used for processing mesostructures comprises: CAD/CAM Titanium base, CEREC®.

CAD/CAM titanium bases, CEREC® are fixated in an implant as prosthetic titanium base for adhesion to mesostructures in order to restore function and aesthetics in the oral cavity.

Contraindications:

- Insufficient oral hygiene
- Insufficient space
- Bruxism
- Single restorations with free-end saddle
- Restorations whose length exceeds the ratio of 1:1.25 compared to the length of the implant
- Restorations that require an angulation correction of more than 15° for the S series and more than 20° for the M and L series to the implant axis.
- We offer tiologic® Ø 3.3 mm implants for patients with narrow alveolar ridges. Due to the small diameter and low load capacity (compared to tiologic® Ø 4.2 mm implants), these implants have a limited range of indications. In fully edentulous cases, four or more tiologic® implants must be inserted with a splinted bar restoration.
- In partially edentulous cases, implant supported restorations must be combined with tiologic® Ø 4.2 mm, Ø 4.8 mm or Ø 5.5 mm implants and a splinted fixed prosthetic restoration.
- In single restorations, tiologic® Ø 3.3 mm implants should only be used for the lower incisors or the upper lateral incisors and only with a minimum 11.0 mm implant length.
- Single restorations on tiologic® Ø 3.7 mm, Ø 4.2 mm, Ø 4.8 mm or Ø 5.5 mm implants require a minimum 9.0 mm implant length.

For the USA only

CAUTION: The federal law of the USA restricts the sale of this product to or on the order of physicians, dentists or licensed practitioners.

5. Notes on processing

5.1 Scanning

- Tighten the CAD/CAM Titanium base, CEREC® on the matching laboratory analog in the master model with the supplied AnoTite screw.
- Place the respective scan body (see section 2) on the CAD/CAM Titanium base, CEREC® making sure it is seamlessly seated using the provided guide groove. The scan body can be scanned without powder or scan spray.
- Scan the model alternatively with inEos Blue, inEos, CEREC 3 or CEREC AC.
- Design the individual shape of the mesostructure with inLab 3D for Abutments V3.65 or with inLab 3D V3.80 (or higher) taking into account the following indications: Since the tiologic® implant system is currently not integrated in the Sirona CEREC® software, the camlog® implant system can be alternatively selected within the software for the design process. Please note: Check the correct position of the restoration before the grinding process. If you choose a camlog® implant system with three positioning options, a correction could be necessary. If this is not observed, the mesostructure could be positioned incorrectly. Correct the position with a corresponding recess in the unsintereed mesostructure moved into the correct position. Then, grind the shape out of an inCoris Zi meso block (see inLab 3D/inLab 3D for Abutments, Operator's Manual). Please observe the indications on designing, processing and sintering zirconium oxide in the instructions for use of inCoris Zi meso blocks.

5.2 Processing the CAD/CAM Titanium base, CEREC

5.2. Processing the CAD/CAM Titanium base, CEREC®
The diameter of the CAD/CAM Titanium base, CEREC® should not be reduced e.g. by grinding. It is not recommended to shorten the CAD/CAM Titanium base, CEREC®. The contact surfaces between the CAD/CAM Titanium base, CEREC® and the implant must not be blasted or processed in any other way.

Only the parts of the CAD/CAM Titanium base, CEREC® surface that will be glued to the mesostructure (50 µm aluminium oxide, max. 2.0 bar) should be blasted and then cleaned with alcohol or steam. Secure the CAD/CAM Titanium base, CEREC® in a laboratory analog to protect the internal connection. Use "PANAVIA™ F 2.0" (adhere to the Instructions for use relating to Panavia: www.kuraray-dental.de) extraoral.

- For better handling during the gluing process, insert the CAD/CAM Titanium base, CEREC® in a laboratory implant or a polishing aid.
- Cover the hexagon socket head of the abutment screw with wax.
- Mix the glue according to the manufacturer's instructions and apply it to the CAD/CAM Titanium base, CEREC®.
- Slide the sintered inCoris Zi mesostructure onto the base as far as it will go. Make sure it latches into the rotation and position stop.
- Remove excess glue immediately.
- Apply the Airblocker ("Oxyguard") to the transition point between ceramic and titanium and to the screw funnel to finish hardening the glue.
- After hardening, remove the residue glue with a rubber polisher.

5.3. Note for dental practitioners

The CAD/CAM Titanium bases are supplied in non-sterile condition. Please observe the instructions of the implant manufacturer.

5.3.1 Sterilization

The individual abutments must be sterilized prior to usage. Please observe the applicable legal regulations in the country and the applicable hygiene standards for the dental practice.

Only use validated sterilization procedures and comply with the hygiene rules for dental practices to sterilize the individual abutments. Observe the sterilization parameters.

Steam sterilization can be performed with the fractional vacuum or the gravitation procedure. The sterilization time is 5 minutes at 132°C/ 270°F or 15 minutes at 121°C/250°F. Perform steam sterilization using only equipment that complies with EN 13060 and EN 285 standards.

The sterilization procedures must be validated in accordance with EN ISO 17664. The responsibility for the sterilization of the individual abutments lies with the user. Make sure to only use suitable equipment, materials and product-specifically validated procedures for sterilization. Use validated procedures only. The equipment must be maintained and serviced regularly and properly. The dental technician that processes the CAD/CAM Titanium base, CEREC® and the mesostructure must inform the dentist of the need to sterilize the components before they are inserted in the patient's mouth!

5.3.2 Tightening torque

AnoTite screw CAD/CAM titanium base on model: manually
in mouth: 30 Ncm

6. Quality / Warranty and Liability / Development

Development, clinical testing, production and quality control of the tiologic® product range are completed in accordance with the Medical Device Directive 93/42/EEC.

Sections 9 and 10 of our General Terms of Delivery and Payment apply with regard to warranty or liability – unless otherwise provided in this quick reference guide.

Warranty and liability are rendered void in particular if the products are not used by the operator or a third party in accordance with the instructions for use; this also applies if the tiologic® product range is used in combination with products of other manufacturers, which have not been specifically recommended for use by Dentaurum Implants.

Dentaurum Implants GmbH has no control over processing and application of the product, which are the sole responsibility of the user.

Technical advice (oral and written) is based on the scientific and technical knowledge available when the product is put on the market. It does not release the user from the responsibility of personally checking the suitability of the products for the intended indication and application. Advice is only given as a non-binding recommendation, which cannot be assumed to provide any form of assurance or guarantee.

All products are subject to continuous development based on current scientific knowledge and we reserve the right to make changes in the construction, design or material of the products.

7. Explanation of symbol

 Please refer to the label. Additional information can be found at www.dentaurum.de (Explanation of the symbols REF 989-313-00).

Date of Information: 05/15

Subject to modifications

CEREC® is a registered trademark of Sirona.

camlog® is a registered trademark of CAMLOG.

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